



TS MANAGER

VALIDATION SYSTEM



Test Support



Sterilization validation




Food Certification



TS Manager is the core of Tecnosoft validation and process monitoring solutions. It is a powerful software that supports all Tecnosoft loggers. It can be used for sterilisation validation and daily routine, mapping, profiling, test support for food certification.

Compatible dataloggers



BlueWave

SterilDisk

MicroW S / L / XL

S-Micro / L / XL

PressureDisk


Pirani Vacuum Logger

RHTemp

TempStick

HumiStick

TC-Log



Create mission profiles, assign each test to customer and instrument such as autoclave, storage, fridge, group data from different loggers, make advanced analysis, print completely customizable reports.



VERSION COMPARISON TABLE

FEATURE	LITE*	BASIC	PRO
Starting & downloading device	●	●	●
Creating profiles for mission programming	●	●	●
Lethality analysis (F0, PU, A0)	●	●	●
Max and Min values analysis	●	●	●
Overkill analysis	●	●	●
Vacuum leak analysis	●	●	●
Give names to devices	●	●	●
View a mission from the archive	●	●	●
Apply analysis to a mission	●	●	●
Mission statistics	●	●	●
Export mission data	●	●	●
Print a mission	●	●	●
Group missions together	●	●	●
Group analysis	●	●	●
Group statistics	●	●	●
Lethality curves on graph	●	●	●
Export data of group	●	●	●
Column's header with logger's name instead of Serial Number	●	●	●
Column's header shortening	●	●	●
Create sections in a group	●	●	●
Create groups from sections	●	●	●
Print group report	●	●	●
Generate new statistics with formula expert	●	●	●
Multi user and multi level	●	●	●
Banned users management	●	●	●
Client/Server installation	●	●	●
Backup and restore function	●	●	●
Customers and instruments database	●	●	●
Audit trail & intergrity check	●	●	●
Data encryption	●	●	●
RHTemp adjustament (on software)	●	●	●
21 CFR Part 11, Annex 11, GAMP compatibility	●	●	●

* Available from Q3 of 2024

TS Manager system (Basic and Pro version only) is compatible with the following standards: ISO 17665, ISO 15883, UNI EN 285